



Low viscosity hydrophilic vinylpolysiloxane for high precision impression in fixed protheses

For professional use in dental clinic

Usage: Recommended for impression as correction fluid

- in the two-stage technique(double impression) in conjunction with a Putty of Extreme range

Characteristic:

- fast polymerization, shorter time in the oral cavity for a greater patient comfort
- high index of flowability
- excellent reproduction of details
- excellent elastic recovery
- high resistance to deformation
- high dimensional stability
- high tear strenght
- hydrophilic - thixotropic
- non-toxic - disinfectable
- violet colour

Packaging : cartridge 50 ml.

Mixing ratio 1:1