



Very high viscosity A-Silicone for fixed protheses.



For professional use in dental clinic

Usage: Recommended as first impression material

- for the two-stage technique (double impression)

Characteristics:

- rapid polymerization, shorter time in the oral cavity for a greater patient comfort
- high index of flowability
- doesn't compress the gingival mucosa
- excellent elastic recovery
- high resistance to deformation
- final hardness 62 Sh. A
- high dimensional stability
- no greasy surface
- light blue colour

Packaging: 2 jars (Catalyst + Base) 300+300 ml.