



Condensation Curing Silicone

Professional use in Dental Laboratory

Professional use in Dental Clinic extra Cee Countries

MICRON PUTTY

MICRON CONFORT

Usage:

- realization of masks in dental Lab.
- matrix support for artificial gengiva
- support for making temporary bridge/crowns
- for taking impression in the two stage technique in association to Micron Light or Micron Regular

Usage:

- realization of masks in dental Lab.
- matrix support for artificial gengiva
- support for making temporary bridge/crowns
- for taking impression in the two stage technique in association to Micron Light or Micron Regular

Characteristics:

- highly millable without crumbling
- excellent flowability
- no deformation
- heat resistance 110°C.
- yellow-green colour
- hardness 62 Sh. A
- mixing with catalyst gel or liquid
- no greasy surface

Characteristics:

- highly millable without crumbling
- excellent flowability
- no deformation
- heat resistance 110°C.
- yellow colour
- hardness 52 Sh. A
- mixing with catalyst gel or liquid
- no greasy surface