



High viscosity addition silicone (polyvinylsiloxane) for a fast and accurate bite registration

For professional use in dental clinic

Usage: Indicated for all intraoral registration, for registration of dental arches and for controls in centric relation.
Compatible with CAD-CAM systems using optical, laser and contact reading.

Characteristic

- Fast polymerization for a major comfort for the patient
- High index of flowability
- Excellent reproduction of the details
- Fast < 60"
- High resistance to deformation
- High dimensional stability
- High rigidity 95 ShA
- Hydrophilic - thixotropic
- Non-toxic - disinfectable
- Colour light blue

Packaging : cartridge 50 ml.

Mixing ratio 1:1